


| | | | |
|----------|---------|------|-----------|
| Customer | Contact | Date | 30-Aug-10 |
|----------|---------|------|-----------|

| | | | | | | | | | |
|---------------------|--|--------------------|-------|----------------------|----------|-------------|----------------------|----------|-------------|
| Date | | Test 2 | | | | | | | |
| Material | Powdered Sugar | Material condition | .5 mm | Bulk Density (gm/cc) | | .1 mm | Bulk Density (gm/cc) | | |
| Desired Sample size | 6 mg tgt. wt. | Patriculate size | Fines | From Table | Measured | Tap Density | From Table | Measured | Tap Density |
| Desired accuracy | Range= 5.4 - +/-10% 6.6 mg (+/- .5 mg) | Crystals >5 micron | >10% | NA | 0.57 | 0.6 | NA | | |

| | | | |
|------------------|------------------------------------|-----------------------------|-----------|
| Desired accuracy | Dispense into 2 ml Wheaton Bottles | rh | 49% |
| Pipette Size | 0.125 dia. | Filter Cup: Filter cup mesh | 10 micron |

| | |
|--------------------|------|
| Micrometer setting | 0.40 |
| Sample weight | Set |
| MG | 5.80 |
| | 5.90 |
| | 5.80 |
| | 5.90 |
| | 5.90 |
| | 6.00 |
| | 6.00 |
| | 5.60 |
| | 5.90 |
| | 6.00 |
| | 5.60 |
| | 6.00 |
| | 6.30 |
| | 6.00 |
| | 6.20 |
| | 6.30 |
| | 6.30 |
| | 5.60 |
| | 5.80 |
| | 6.10 |
| | 6.30 |
| | 5.70 |
| | 5.80 |

| | | |
|-----------------------|----------------|-----------|
| Control unit settings | Vacuum (in hg) | Air (psi) |
| | 12 | 14 |




| | | | | |
|------------------|----------|-------|----------|-------|
| Cycle time (sec) | Aspirate | Level | Dispense | total |
| | 1 | 1 | 1 | 3 |

| | |
|--------|---------|
| 5.95 | Av |
| 5.95 | Mean |
| 0.2337 | Std Dev |
| 6.3 | High |
| 5.8 | Low |
| 7.94 | % |

Range

| | |
|---------|---------|
| #DIV/0! | Av |
| #NUM! | Mean |
| #DIV/0! | Std Dev |
| | High |
| | Low |
| | % |

Range



Notes/observations

Powdered sugar (as an analog to the sodium cyanoborohydrite)
 1a. White crystals, easily aspirated and dispensed. No issues dispensing into the customers bottles.
 1b. Set up and calibration are straight forward.
 1c. Requires relatively high air pressure to eject sample from the pipette due to the crystalline shape compacting.
 1d. Tip pressure on leveling plate compacts the crystals so the sugar is expelled as a loose pellet.
 Care must be exercised to ensure all powder proud of the tip orifice must be removed to minimize sample weight variation. See samples # 11 - 15 and 16 - 21.
 1f. No special process efforts were taken to optimize the dispense weight variation as this powder is an analog and may vary in characteristics from the actual powder.
 2.0 Standard 0.125 dia tip used for testing due to small dispense volume.
 3.0 **Aerosol:** No aerosol or loose crystals were observed during or after dispense.
 4.0 Under normal laboratory conditions, periodic cleaning, based on humidity. Tip/filter cleaning and changeover is <5 minutes.

Summary: Easy material to pipette, even at the small dispense weight. This powder is an analog and may vary in characteristics from the actual powder.

