

Customer							
Material	Moly Tri-Oxide MoO3		Est. Bulk Density	0.5			
Desired range	.15 gm plus/minus 0.0002 gm						
Desired accuracy	Plus/Minus 0.0002 gm						
		Sample weight	Control unit settings				
Micrometer setting		gm	Vacuum (1in hg)	Air (psi)			
500	Start	0.1252	10	10	low pack		
		0.1252					
600		0.1684					
550	Run	0.1461	22	17	Reset Air/Vac		
		0.1592					
		0.1505					
		0.1596					
		0.1550					
		0.1564					
		0.1569					
		0.1594					
		0.1618					
		0.1624					
		0.1520					
		0.1526					
		0.1540					
		0.1558 Av					
		0.1558 Mean					
		0.0048 STD DEV					

