

Customer						
Material	PEANUT NUT MEAT	Est. Bulk Density	0.5	Date	4/23/02	
Desired Sample Size	100 mg	Pipette Size	.250 dia	rh	50%	
Desired accuracy	Plus/Minus 5 mg					
Micrometer setting	Sample weight gm	Control unit settings Vacuum (in hg) Air (psi)				
2.75	106.8	7	10	Material is damp		
	106.7			↓		
	105.1			Sample drying		
	109.2			↓		
	108.4					
	103.5					
	102.8					
	100.1					
	99.7					
	99.5					
	98.5					
	98.1					
	103.2 Av					
	103.1 Mean					
	3.997 Std Dev					
Material tends to clump. Oily consistency						

