



PIPETTER VOLUME RAW DATA		PIPETTE S/N	SAMPLE MATERIAL										DATE	SAMPLE MATERIAL = TALC	
PIPETTER ORICE DIAMETER	SAMPLE NUMBER	MICROMETER THIMBLE SETTING											Rh	?	
.500 dia.	gm	1.0000	2	3	4	5	6	7	8	9	10	1.3	VAC / AIR	COMMENTS	
	1	0.6847				1.8776					3.2973	4.1280	5 in. hg /6 PSI		
	2	0.6842				1.8992					3.2610	4.1440			
	3	0.6846				1.8996					3.2978	4.1364			
	4	0.6833				1.8867					3.2579	4.1930			
	5	0.6830				1.9203					3.2675	4.1815			
	6	0.6850									3.2892	4.2225			
	7	0.6840				1.8967	Av				3.2543	4.1648			
	8	0.6835				1.8966	Mean					4.2145			
	9	0.6855				0.0161	Std. Dev.			Av	3.2750	4.1258			
	10	0.6852								Mean	3.2749	4.1530			
	11	0.6860								Std. Dev.	0.0191	4.1807			
	12														
	13	0.6845	Av									4.1677	Av		
	14	0.6845	Mean									4.1675	Mean		
	15	0.0010	Std. Dev.									0.0336	Std. Dev.		
	16														



27														
28														
29														
30														
31														